

VES MATIC 5

Italian Design meets artificial intelligence



US FDA
Registered



**Continuous and
Random Access
Analyzer.**



**Maximum
Operator safety
and Protection.**



**Higher throughput
120-190 test
results/hour**



**Modern and
Interactive Graphic
User Interface.**

Key Benefits



Operational

Dry Technology

- User Friendly touch screen
- First result is obtained in 28 minutes, subsequent every sample is reported in 19 seconds.
- Standardization in mixing of samples with Inbuilt mixer
- Inbuilt printer for documentation of report
- Compatible with all types of racks and all types of tubes
- Archival of upto 25000 sample results



Clinicians

- High co-relation with reference 'Gold Standard' (ICSH Approved)
- Automated HCT correction facility eliminates manual calculation for Anaemic sample.
- Artificial intelligence based recognition of lipemic, hemolyzed and coagulated samples.



Patient

- No dedicated tubes required for sample analysis. Less sample volume required



Economical

- No reagent and calibrator



Society

- No waste generation

Ordering Information

| Product Code | Description |
|--------------|------------------------------------|
| 166823 | VES MATIC 5 ESR Analyzer |
| 165941 | Next 5 K Transponder (5000 Test) |
| 165942 | Next 10 K Transponder (10000 Test) |
| 195086 | ESR Control Cube (2X9 ml) |

Marketed by

Transasia Bio-Medicals Ltd.

An ISO 13485:2016 / 9001:2015 Certified Company

Transasia House, 8, Chandivali Studio Road,
Chandivali, Andheri (E), Mumbai - 400 072

T: (022) 4030 9000 F: (022) 2857 3030

✉ responses@transasia.co.in