

Erba LisaWash[®] - II

Automated ELISA Microplate Washer



50+

50 user defined protocols

<2 μ L

Residual volume of <2 μ L/ well



Automatic monitoring of vacuum pressure to avoid clogs or blockages

Making Automation Affordable
for labs Everywhere

LisaWash II ELISA Plate washer was created to address the need for effective washing, a fundamental step to any successful ELISA. It is a Fully automated, easy operating system with a user friendly software. Equipped with an easily removable 8 or 12 way manifolds Erba LisaWash is compatible for all types of microplates or microstrip and has a High-efficiency ELISA washing with automated monitoring of Vacuum. This enhances the accuracy of washing drastically providing assurance to the users and provides effective washing for all kinds of ELISA assays.

Key Benefits



Dependability

- LisaWashII is a reliable, simple to use ELISA Washer designed for automated and reproducible/efficient washing of 96 well plates.
- 2 point aspiration per well ensures minimal residual volume, which eliminates risk of high background signal to ensure reliable results without compromising on the analytical sensitivity.



Flexibility

- Washing performance for various plate geometries can be optimized by adjusting various wash parameters like, cycles, soak time, Top & bottom well adjustments, dispense pressure, provides secure washing performance without risk of overflow, for routine as well as research application.



Simplicity

- LisaWashII is so simple to set up that anyone can do it. Easy setup of wash protocols, aided by a simple intuitive user interface requiring minimal training, makes it very convenient to use.



Safety

- Waste sensor helps ensure safe operations & user safety, apart from preventing backflow & damage to electronics. Simple structure ensures that decontamination can be carried out effortlessly without engineering support

Ordering Information

Product Code	Description
165530	Lisa Wash II

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno, Czech Republic

Phone: +420 517 077 111

✉ sales@erbamannheim.com

www.erbamannheim.com

Transasia Bio-Medicals Ltd.

An ISO 13485:2016 / 9001:2015 Certified Company

Transasia House, 8, Chandivali Studio Road,

Chandivali, Andheri (E), Mumbai - 400 072

T: (022) 4030 9000 F: (022) 2857 3030

✉ responses@transasia.co.in

