

ECL 412

Four Channel Coagulation Analyzer
for Small to Mid size labs



High
Performance



Reliable and
Accurate Results



User-friendly
software

Making Automation Affordable
for labs Everywhere

Transasia's ECL 412 is semi-automated Coagulation analyzer with four testing position providing operators with a compact easy to use system. ECL 412 is a blend of modernity & simplicity. Transasia also offers a reagent portfolio as a complete way to your preference in Hemostasis.

Key Benefits



Operational

- High Performance: Test performed are PT, APTT, FIB, Intrinsic and Extrinsic Factor Assay, LA, PC, PS, AT III, D-dimer
- Reliable and Accurate result- Based on gold standard Nephelometry Optics Principle
- User-friendly software - Large graphic alpha numeric color touch screen display
- Complete QC programming- Multilevel Controls with L-J graph
- Reduced Waiting Time- Built-in incubator with 20 reaction cuvette position
- Automated Results- Microprocessor based analyzer
- Audit trail for Reagents, Controls, Patients sample results
- Elimination of transcription errors through LIS connectivity



Clinicians

- Helps in diagnosis of Bleeding and Thrombosis disease
- Accurate results and processing of samples with endogenous interference as Hemolysis, Icteric and Lipemic (HIL)
- Built-in Graphical Thermal Printer



Economical

- Labs can offer coagulation Testing for Pre and Post Surgery, Cardiac, Gynec patients

Erba Coagulation Line

Instruments	Reagents	
ECL 105	Erba Protime LS	Erba D Dimer
ECL 412	Erba Actime	Erba Standard Plasma
ECL 760	Erba Calcium Chloride	Erba control N/P
	Erba Thrombin Reagent	Erba control N/P Plus

Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno, Czech Republic

Phone: +420 517 077 111

✉ sales@erbamannheim.com

www.erbamannheim.com

Transasia Bio-Medicals Ltd.

An ISO 13485:2016 / 9001:2015 Certified Company

Transasia House, 8, Chandivali Studio Road,

Chandivali, Andheri (E), Mumbai - 400 072

T: (022) 4030 9000 F: (022) 2857 3030

✉ responses@transasia.co.in

