


New

ErbaQik[®]

Rapid Tests for Dengue & Malaria

with

DUAL | COLOUR
ADVANTAGE



Simple, Accurate & Convenient!

ErbaQik[®]




Rapid Tests for Dengue & Malaria


Innovation

DUAL COLOUR ADVANTAGE

ErbaQik rapid test utilizes unique **Black Gold Particle**, which makes it the only Bi-colour detection Immunochromatographic device available for the diagnosis of Malaria/Dengue infection.

Benefits:

-  Ease of interpretation
-  Reliable results
-  Unmatched sensitivity and specificity - Best in Class*



Black Gold Particle

ErbaQik[®] Malaria Pan Ag

(Pan pLDH)

Sample: Whole Blood

Sensitivity: 97.14%

Specificity: 100%

ErbaQik[®] Malaria Pf/Pan Ag

(HRP2/Pan pLDH)

Sample: Whole Blood

Sensitivity: 97.14%

Specificity: 100%

ErbaQik[®] Malaria Pf/Pv Ag

(HRP2/Pv pLDH)

Sample: Whole Blood

Sensitivity: 97.14%

Specificity: 100%

ErbaQik[®] Dengue NS1 Ag

(Dengue NS1 Ag)

Sample: Serum/Plasma/WB

Sensitivity: 93.75%

Specificity: 100%

ErbaQik[®] Dengue IgG/IgM

(Dengue IgG/IgM Antibodies)

Sample: Serum/Plasma/WB

Sensitivity: IgG-100%, IgM-92.59%

Specificity: IgG & IgM- 100%

ErbaQik[®] Dengue Duo

(NS1 Ag + IgG/IgM Ab)

Sample: Serum/Plasma/WB

Sensitivity: NS1-93.75%, IgG-100%, IgM-92.59%

Specificity: NS1-100%, IgG-100%, IgM-100%

To place orders:

 1800 103 8226

 91369 50470

* As per the external evaluations done at NABL labs



ERBA Diagnostics Mannheim GmbH

Mallastrasse 69-73 68219 Mannheim, Germany
Telephone : (+49) 621 8799770 Fax : (+49) 621 8799688
Email : sales@erbamannheim.com

TRANSASIA BIO - MEDICALS LTD.

Transasia House, 8 Chandivali Studio Road, Andheri (E), Mumbai - 400 072
Tel.: (022) 4030 9000, Fax: (022) 2857 3030
Email: responses@transasia.co.in

Stay connected:    

www.erbamannheim.com