

LAURA XL

Fully Automated Urine Chemistry
and Sediment Analyzer



Faster turn- around time

Achieve quickest
results every time



Reduced manual work

Optimise your
work-flow



Image storage

See & Share more
with ease

Making Automation Affordable
for labs Everywhere

Laura XL is complete walkaway system combining urine strip reading and digital microscopy on a single platform while making it a preferred choice for high workload laboratories.

Key Benefits



Operational

- Manual Work Minimization: The fully automated aspiration, processing & analysis of Urine sample would reduce required manual work in slide preparation & viewing under microscope
- Error Reduction: The objective strip reading, software based sediment analysis & automated result transmission could possibly minimize the amount of errors with manual work
- Time Saving: With the complete walkway system it could help use valuable technician time for other important work
- Reduced Review Rate: The On-screen images of strip & sediment could be referred for any repeat analysis if required
- Throughput of 125 test/hr in Sediment/Hybrid mode & 160 test/hr in Chemistry mode
- Memory of 5 lakh results
- Sample volume required of 2 ml urine
- Sample loading capacity of 120
- Internal barcode scanning facility & LIS compatibility



Clinical

- Advanced Reporting: The facility to save & share high resolution sediment images shall be helpful to prove the analysis
- Auto detection of Color & Clarity
- Chemistry parameters include: Glucose, Protein, pH, Specific Gravity, Leucocytes, Blood, Bilirubin, Urobilinogen, Ketone & Nitrite (DekaPHAN AUTO)
- Sediment parameters include: WBC, WBC clumps, RBC (with dysmorphic sub-classification), Epithelial cells (Squamous & Non-Squamous), Casts (Hyaline & Pathological), Crystals (Uric acid, Triple Phosphate, Calcium oxalate & Small crystals), Bacteria (Cocci & Rods), Yeast, Mucous, Sperm

Ordering Information

Product Code	Product	Quantity
165467	LAURA XL	1 unit
132419	DekaPHAN AUTO	100 strips
132399	OptiSOL XL 1500	400 ml
133506	OptiSOL 750	6 x 50 ml
134450	Urinorm XL	2 x 120 ml

Erba Lachema s.r.o.

Karásek 2219/Id, 621 00 Brno, Czech Republic
Phone: +420 517 077 111

✉ sales@erbamannheim.com

www.erbamannheim.com

Transasia Bio-Medicals Ltd.

An ISO 13485:2016 / 9001:2015 Certified Company

Transasia House, 8, Chandivali Studio Road,
Chandivali, Andheri (E), Mumbai - 400 072

T: (022) 4030 9000 F: (022) 2857 3030

✉ responses@transasia.co.in

