

# XL-1000

Fully Automated Random Access Clinical Chemistry Analyzer  
*A perfect solution for high workload Laboratories*



**US FDA  
REGISTERED**



**Increase  
Productivity**

Maximise testing  
capabilities to next  
level



**Accurate  
Results**

Achieve highest  
accuracy first time,  
every time



**Quality  
First**

QC capabilities  
put patient safety  
first

Making Automation Affordable  
for labs Everywhere

XL-1000 is a random access, discrete, patient prioritized, fully automated clinical chemistry analyzer with autoloader designed to meet the needs of today's modern clinical laboratories with high workload. Once samples are loaded, it is a complete walk-away system.

## Key Benefits



### Flexibility and versatility

- Throughput: 800 tests / hour photometric and 1040 tests / hour with ISE (optional Direct ISE with Na+, K+, Cl-, Li+)
- 86 positions for reagents with on-board cooling
- 159 sample positions (including blank, controls, calibrators, upto 44 stat samples) with sample autoloader rack system (continuous loading sample racks)



### Quality first

- QC capabilities put patient safety at the centre of your laboratory
- 2 level QC plot data with QC rules, provision for lab mean and Twin plot
- QC for Urine parameters



### Economical in all aspects

- 147 permanent glass cuvettes helps to reduce CPRT
- Reaction volume only 150 µl ensures less reagent consumption for every test
- < 25 litres water consumption



### High accuracy and precision

- Separate probes for R1 and R2 reagents and sample with internal and external washing ensures no carry over
- Dual mixer with variable speed technology ensures proper mixing of reagent and samples
- Multi-wavelength diffraction grating with 15 wavelengths (340, 376, 405, 450, 480, 505, 528, 546, 570, 600, 628, 660, 700, 750, 800 nm) ensures precise reading every time.
- Clot detection function, and serum Indices identification provides next level sample check
- 10 stage Twin onboard laundry ensure no cross contamination



### Smart functions

- Intelligent remote service support using IoT for reporting errors for quick service
- Inbuilt barcode for reagent and samples
- User defined reflex testing
- Touch screen with all new user-friendly GUI
- Bi-Directional LIS

## Ordering Information

Product Code	Product	Description
190060	XL 1000	FAA - Floor Type

### Erba Lachema s.r.o.

Karásek 2219/1d, 621 00 Brno, Czech Republic

Phone: +420 517 077 111

✉ sales@erbamannheim.com

[www.erbamannheim.com](http://www.erbamannheim.com)

### Transasia Bio-Medicals Ltd.

**An ISO 13485:2016 / 9001:2015 Certified Company**

Transasia House, 8, Chandivali Studio Road,

Chandivali, Andheri (E), Mumbai - 400 072

T: (022) 4030 9000 F: (022) 2857 3030

✉ responses@transasia.co.in

