


E r b a **Lisa**<sup>®</sup>  
Thyroid Hormone ELISAs



Quantitative **ELISA** for Thyroid Hormones


Erba **Lisa**® Thyroid ELISA kits are methodically designed for the accurate quantitative measurement of thyroid hormones in human serum.

Erba **Lisa**®  
Total Triiodothyronine  
(tT3)




Sensitivity: 0.066 ng/mL

Erba **Lisa**®  
Total Thyroxine  
(tT4)



Sensitivity: 0.32 µg/dL

Erba **Lisa**®  
Thyroid Stimulating Hormone  
(TSH)



Sensitivity: 0.045 µIU/mL



Product Code	Product Name	Pack Size
133321	Erba <b>Lisa</b> ® Total Triiodothyronine (tT3)	96 Test
133322	Erba <b>Lisa</b> ® Total Thyroxine (tT4)	96 Test
132904	Erba <b>Lisa</b> ® Thyroid Stimulating Hormone (TSH)	96 Test

### Compatible ELISA Instruments:

**ELAN 30s**



**Erba LisaScan® EM**



**Erba LisaWash® - II**



**ERBA Diagnostics Mannheim GmbH**

Mallaustrasse 69-73 68219 Mannheim, Germany  
Telephone : (+49) 621 8799770 Fax : (+49) 621 8799688  
Email : sales@erbamannheim.com

Toll Free Number

**1800 103 8226**

W: www.erbamannheim.com  
Stay connected: [in](#) [t](#) [e](#) [f](#)

**TRANSASIA BIO-MEDICALS LTD.**

Transasia House, 8 Chandivali Studio Road,  
Andheri (E), Mumbai - 400 072. T: (022) 4030 9000,  
F: (022) 2857 3030 E: responses@transasia.co.in

