

Erba Chem 5[®]x

- Extensive Biochemistry test menu
- Elaborate Q.C. Program and Memory
- Inbuilt Printer
- Graphic Display
- 5000 Test Result storage



Clinical Chemistry Analyzer



Erba Chem 5[®]x

Clinical Chemistry Analyzer

TECHNICAL SPECIFICATIONS

- 150 open user-defined test programs
- 9 Assay modes : 1 point Linear, 2 point Linear, 1 point Non Linear, 2 point Non-linear, Rate A Linear, Rate A Non-Linear, 1 point Sample Blank Linear, 1-point Sample Blank Non-Linear, Absorbance
- Mono & Bichromatic measurement
- Photometric range from 0 to 3.0 OD
- Photometric resolution up to 0.0001
- Automatic Zero Setting
- Static Photometer (340-630 nm)
- 6 Wavelengths: 340, 405, 505, 546, 578, 630 nm + 2 Optional Positions
- Multi-Point calibration with graphical display
- Levey-Jennings graph for quality control
- 5000 Test Patient Results storage
- Peristaltic Pump Aspiration System with calibration facility
- Dual Cuvette System : Flow cell volume - 33 µl
: 10 mm external glass / Polystyrene
- Aspiration Volume from 350 to 999 µl
- Temperature Control at 25°, 30° & 37°
- Repeat measurement by Read key
- Collated report date wise, ID wise (Alphanumeric)
- Inbuilt Thermal Printer
- Direct Printer Connectivity through USB Port
- USB Connectivity to PC
- Connectivity to external Incubator (Optional with 16 position)

Specifications are subject to change without prior notification.

Distributor Stamp

ERBA Diagnostics Mannheim GmbH

Mallaustrasse 69-73 68219 Mannheim, Germany
Telephone : (+49) 621 8799770 Fax : (+49) 621 8799688
Email : sales@erbamannheim.com

TRANSASIA BIO-MEDICALS LTD.

Transasia House, 8 Chandivali Studio Road, Andheri (E), Mumbai - 400 072
Tel.: (022) 4030 9000, Fax : (022) 2857 3030
Email : responses@transasia.co.in

Stay connected:     

www.erbamannheim.com

